



Detergent tank heater

Unit for precise heating up to max. 80 °C of cleaning media in the cleaning tank.



Part numbers:

- 090520-6HCT (AF20)
- 090530-6HCT (AF30)
- 0905AQ-6HCT (AQUBE® / PSE 5, 6, 7, MV3 ONE)
- 0905AQ8-6HCT (AQUBE® / PSE MV8)
- 0905AQ9-6HCT (AQUBE® / PSE 9)

Application

Certain detergents develop their optimal cleaning effect only at a certain temperature. When using such media, it is important to ensure their optimum temperature precisely in order to obtain reproducibly perfect cleaning results.

The cleaning temperature of the unit equipped with a contact protection will be preset to a defined value. The PLC control module for the heater enables the activation of the heating element, monitors the temperature and activates the built-in overtemperature protection.

kolb detergents in general are high-performance low-temperature cleaners, which generally reach their optimum efficiency already at room temperature (about 20 °C) or when heated to a maximum of 45 °C.

Compatibility

AQUBE® LH5
AQUBE® LH6
AQUBE® LH7
AQUBE® XH6
AQUBE® XH7
AQUBE® XH7 TZD
AQUBE® XH9
AQUBE® XV7
AQUBE® XV9
AQUBE® MH6
AQUBE® MH7
AQUBE® MH7 TZD
AQUBE® MH9
AQUBE® MV3 ONE
AQUBE® MV8 sTWIN
AQUBE® MV9 QUAD
PSE LH5
PSE LH6
PSE LH7
PSE MH6
PSE MH7
PSE MV8 TWIN
PSE 300 MV3 XL
PSE 300 MV3 XXL
PSE 300 MV3 ONE
PSE 300 VL
PSE 300 VS
PSE 300 HD
AF 04
AF 20
AF30
AF75

Parts list

- 1x screw heater device / Angulated heating rod
- SPS heating control add on module
- 1x contact protection
- 1x blending unit (only with AF20 und AF30 systems)

Standard in the system

Option ex factory

Retrofittable by customer

Not compatible